

## Datasheet

### c-Myc tag monoclonal antibody, clone 3E8

**Catalog Number:** MAB15241

**Regulatory Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against synthetic peptide of human c-Myc tag.

**Clone Name:** 3E8

**Immunogen:** A synthetic peptide corresponding to full length human c-Myc tag.

**Host:** Mouse

**Applications:** IF, IP, WB-Re  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** The antibody can detect c-Myc fusion proteins.

**Form:** Liquid

**Purification:** Affinity purification

**Isotype:** IgG1

**Recommend Usage:** Immunofluorescence (1:1000)  
Immunoprecipitation (1:200)  
Western Blot (1:5000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS (50% glycerol, 0.5% BSA, 0.02% sodium azide).

**Storage Instruction:** Store at -20°C.  
Aliquot to avoid repeated freezing and thawing.