

Datasheet

TFRC monoclonal antibody, clone FG2/12 (PerCP)

Gene Alias: CD71, TFR, TFR1, TRFR

Catalog Number: MAB15288

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against human TFRC.

Clone Name: FG2/12

Immunogen: Human TFRC.

Host: Mouse

Theoretical MW (kDa): 91

Reactivity: Human

Applications: Flow Cyt

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PerCP

Purification: Affinity chromatography

Isotype: IgG1

Recommend Usage: Flow Cytometry (20 μ L/ 10^6 cells)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).

Storage Instruction: Store at 4°C. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 7037

Gene Symbol: TFRC