

Datasheet

Cd19 monoclonal antibody, clone 1D3 (APC)

Catalog Number: MAB15294

Regulatory Status: For research use only (RUO)

Product Description: Rat monoclonal antibody raised
against mouse Cd19.

Clone Name: 1D3

Immunogen: Mouse Cd19-transfected cell line.

Host: Rat

Reactivity: Mouse

Applications: Flow Cyt
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Form: Liquid

Conjugation: APC

Purification: Affinity chromatography

Isotype: IgG2a