

Datasheet

Cd19 monoclonal antibody, clone 1D3 (PerCP-Cyanine5.5)

Catalog Number: MAB15297

Regulatory Status: For research use only (RUO)

Product Description: Rat monoclonal antibody raised against mouse Cd19.

Clone Name: 1D3

Immunogen: Mouse Cd19-transfected cell line.

Host: Rat

Reactivity: Mouse

Applications: Flow Cyt
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PerCP-Cyanine5.5

Purification: Affinity chromatography

Isotype: IgG2a