

Datasheet

Cd8a monoclonal antibody, clone 53-6.72 (PE)

Catalog Number: MAB15301

Regulatory Status: For research use only (RUO)

Product Description: Rat monoclonal antibody raised against mouse Cd8a.

Clone Name: 53-6.72

Immunogen: Mouse thymocytes/splenocytes.

Host: Rat

Reactivity: Mouse

Applications: Flow Cyt
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PE

Purification: Affinity chromatography

Isotype: IgG2b