

## Datasheet

### **Cd5 monoclonal antibody, clone 53-7.3 (PE)**

**Catalog Number:** MAB15308

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rat monoclonal antibody raised  
against mouse Cd5.

**Clone Name:** 53-7.3

**Immunogen:** Mouse thymocytes/splenocytes.

**Host:** Rat

**Reactivity:** Mouse

**Applications:** Flow Cyt  
(See our web site product page for detailed applications  
information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product  
page for detailed protocols

**Form:** Liquid

**Conjugation:** PE

**Purification:** Affinity chromatography

**Isotype:** IgG2a