

Datasheet

Cd5 monoclonal antibody, clone 53-7.3 (PerCP)

Catalog Number: MAB15309

Regulatory Status: For research use only (RUO)

Product Description: Rat monoclonal antibody raised
against mouse Cd5.

Clone Name: 53-7.3

Immunogen: Mouse thymocytes/splenocytes.

Host: Rat

Reactivity: Mouse

Applications: Flow Cyt
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Form: Liquid

Conjugation: PerCP

Purification: Affinity chromatography

Isotype: IgG2a