

Datasheet

Il7r monoclonal antibody, clone A7R34 (FITC)

Catalog Number: MAB15315

Regulatory Status: For research use only (RUO)

Product Description: Rat monoclonal antibody raised against mouse Il7r.

Clone Name: A7R34

Immunogen: Il7r-IgG1 fusion protein.

Host: Rat

Reactivity: Mouse

Applications: Flow Cyt

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: FITC

Purification: Affinity chromatography

Isotype: IgG2a