

## Datasheet

### **Cd4 monoclonal antibody, clone GK1.5 (PerCP)**

**Catalog Number:** MAB15322

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rat monoclonal antibody raised against mouse Cd4.

**Clone Name:** GK1.5

**Immunogen:** Mouse CTL clone V4.

**Host:** Rat

**Reactivity:** Mouse

**Applications:** Flow Cyt  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** PerCP

**Purification:** Affinity chromatography

**Isotype:** IgG2b