

## Datasheet

### **Cd44 monoclonal antibody, clone IM7 (FITC)**

**Catalog Number:** MAB15325

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rat monoclonal antibody raised against mouse Cd44.

**Clone Name:** IM7

**Immunogen:** Dexamethasone-induced cells of the SJL mouse spontaneous myeloid leukemia M1 .

**Host:** Rat

**Reactivity:** Mouse

**Applications:** Flow Cyt

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** FITC

**Purification:** Affinity chromatography

**Isotype:** IgG2b