

Datasheet

Cd44 monoclonal antibody, clone IM7 (PerCP)

Catalog Number: MAB15327

Regulatory Status: For research use only (RUO)

Product Description: Rat monoclonal antibody raised against mouse Cd44.

Clone Name: IM7

Immunogen: Dexamethasone-induced cells of the SJL mouse spontaneous myeloid leukemia M1 .

Host: Rat

Reactivity: Mouse

Applications: Flow Cyt

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PerCP

Purification: Affinity chromatography

Isotype: IgG2b