

## Datasheet

### Il2ra monoclonal antibody, clone PC61 (FITC)

**Catalog Number:** MAB15333

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rat monoclonal antibody raised against mouse Il2ra.

**Clone Name:** PC61

**Immunogen:** IL2-dependent cytolytic mouse T-cell clone B6.1.

**Host:** Rat

**Reactivity:** Mouse

**Applications:** Flow Cyt  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** FITC

**Purification:** Affinity chromatography

**Isotype:** IgG1