

Datasheet

Il2ra monoclonal antibody, clone PC61 (PerCP)

Catalog Number: MAB15335

Regulatory Status: For research use only (RUO)

Product Description: Rat monoclonal antibody raised against mouse Il2ra.

Clone Name: PC61

Immunogen: IL2-dependent cytolytic mouse T-cell clone B6.1.

Host: Rat

Reactivity: Mouse

Applications: Flow Cyt
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PerCP

Purification: Affinity chromatography

Isotype: IgG1