

Datasheet

SELL monoclonal antibody, clone HI62L (PerCP)

Catalog Number: MAB15355

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against human SELL.

Clone Name: HI62L

Immunogen: Peripheral blood mononuclear cells.

Host: Mouse

Theoretical MW (kDa): 76

Reactivity: Human

Applications: Flow Cyt
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PerCP

Purification: Affinity purification

Isotype: IgG2a

Recommend Usage: Flow Cytometry (20 μ L/ 10^6 cells)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).

Storage Instruction: Store in the dark at 4°C.

Entrez GeneID: 6402

Gene Symbol: SELL

Gene Alias: CD62L, LAM-1, LAM1, LECAM1, LNHR, LSEL, LYAM1, Leu-8, Lyam-1, PLNHR, TQ1, hLHRc

Gene Summary: SELL is a cell surface component that is a member of a family of adhesion/homing receptors which play important roles in leukocyte-endothelial cell interactions. The molecule is composed of multiple domains: one homologous to lectins, one to epidermal growth factor, and two to the consensus repeat units found in C3/C4 binding proteins. [provided by RefSeq]