

Datasheet

HLA-DR monoclonal antibody, clone GRB-1 (PerCP)

Catalog Number: MAB15425

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against native human HLA-DR.

Clone Name: GRB-1

Immunogen: Mononuclear cell leukemia acute undifferentiated.

Host: Mouse

Theoretical MW (kDa): 36

Reactivity: Human

Applications: Flow Cyt
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PerCP

Purification: Affinity purification

Isotype: IgG2a

Recommend Usage: Flow Cytometry (20 μ L/ 10^6 cells)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).

Storage Instruction: Store in the dark at 4°C. Avoid prolonged exposure to light.