

## Datasheet

### **Cd28 monoclonal antibody, clone 37.51 (FITC)**

**Catalog Number:** MAB15432

**Regulation Status:** For research use only (RUO)

**Product Description:** Hamster monoclonal antibody raised against native mouse Cd28.

**Clone Name:** 37.51

**Immunogen:** Mouse EL-4 (T-cell lymphoma) cells.

**Host:** Hamster

**Theoretical MW (kDa):** 44

**Reactivity:** Mouse

**Applications:** Flow Cyt  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** FITC

**Purification:** Affinity purification

**Isotype:** IgG2

**Recommend Usage:** Flow Cytometry  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).

**Storage Instruction:** Store in the dark at 4°C. Avoid prolonged exposure to light.

**Entrez GeneID:** 12487

**Gene Symbol:** Cd28

**Gene Alias:** -