

Datasheet

Cd28 monoclonal antibody, clone 37.51 (PerCP)

Catalog Number: MAB15434

Regulation Status: For research use only (RUO)

Product Description: Hamster monoclonal antibody raised against native mouse Cd28.

Clone Name: 37.51

Immunogen: Mouse EL-4 (T-cell lymphoma) cells.

Host: Hamster

Theoretical MW (kDa): 44

Reactivity: Mouse

Applications: Flow Cyt

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PerCP

Purification: Affinity purification

Isotype: IgG2

Recommend Usage: Flow Cytometry

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).

Storage Instruction: Store in the dark at 4°C. Avoid prolonged exposure to light.

Entrez GeneID: 12487

Gene Symbol: Cd28

Gene Alias: -