

## Datasheet

### **Cd45.2 monoclonal antibody, clone 104 (CF-Blue)**

**Catalog Number:** MAB15444

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against native mouse Cd45.2.

**Clone Name:** 104

**Immunogen:** B10.S mouse thymocytes and splenocytes.

**Host:** Mouse

**Theoretical MW (kDa):** 180-240

**Reactivity:** Mouse

**Applications:** Flow Cyt  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** CF-Blue

**Purification:** Affinity purification

**Isotype:** IgG2a

**Recommend Usage:** Flow Cytometry  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).

**Storage Instruction:** Store in the dark at 4°C. Avoid prolonged exposure to light.