

Datasheet

Cd45.2 monoclonal antibody, clone 104 (PerCP-Cyanine5.5)

Catalog Number: MAB15445

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against native mouse Cd45.2.

Clone Name: 104

Immunogen: B10.S mouse thymocytes and splenocytes.

Host: Mouse

Theoretical MW (kDa): 180-240

Reactivity: Mouse

Applications: Flow Cyt
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PerCP-Cyanine5.5

Purification: Affinity purification

Isotype: IgG2a

Recommend Usage: Flow Cytometry
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).

Storage Instruction: Store in the dark at 4°C. Avoid prolonged exposure to light.