

## Datasheet

### H2-D1 monoclonal antibody, clone 34-2-12 (PerCP)

Gene Alias: H-2D, H2-D

**Catalog Number:** MAB15449

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against native mouse H2-D1.

**Clone Name:** 34-2-12

**Immunogen:** (C57BL/6 x DBA/2) F1 mouse splenocytes.

**Host:** Mouse

**Theoretical MW (kDa):** 40

**Reactivity:** Mouse

**Applications:** Flow Cyt  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** PerCP

**Purification:** Affinity purification

**Isotype:** IgG2a

**Recommend Usage:** Flow Cytometry  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).

**Storage Instruction:** Store in the dark at 4°C. Avoid prolonged exposure to light.

**Entrez GeneID:** 14964

**Gene Symbol:** H2-D1