

Datasheet

Cd3e monoclonal antibody, clone 145-2C11 (PerCP)

Gene Alias: AI504783, CD3, CD3epsilon, T3e

Catalog Number: MAB15454

Regulation Status: For research use only (RUO)

Product Description: Hamster monoclonal antibody raised against native mouse Cd3e.

Clone Name: 145-2C11

Immunogen: H-2Kb specific cytotoxic T lymphocyte clone BM10-37.

Host: Hamster

Theoretical MW (kDa): 25

Reactivity: Mouse

Applications: Flow Cyt
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PerCP

Purification: Affinity purification

Isotype: IgG

Recommend Usage: Flow Cytometry
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).

Storage Instruction: Store in the dark at 4°C. Avoid prolonged exposure to light.

Entrez GeneID: 12501

Gene Symbol: Cd3e