

## Datasheet

### **Cd80 monoclonal antibody, clone 16-10A1 (FITC)**

TS/A-1

**Catalog Number:** MAB15457**Regulation Status:** For research use only (RUO)**Product Description:** Hamster monoclonal antibody raised against native mouse Cd80.**Clone Name:** 16-10A1**Immunogen:** Mouse Cd80 (B7) transfected cell line.**Host:** Hamster**Theoretical MW (kDa):** 60**Reactivity:** Mouse**Applications:** Flow Cyt

(See our web site product page for detailed applications information)

**Protocols:** See our web site at<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols**Form:** Liquid**Conjugation:** FITC**Purification:** Affinity purification**Isotype:** IgG2**Recommend Usage:** Flow Cytometry

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).**Storage Instruction:** Store in the dark at 4°C. Avoid prolonged exposure to light.**Entrez GeneID:** 12519**Gene Symbol:** Cd80**Gene Alias:** B7-1, B7.1, Cd28l, Ly-53, Ly53, MIC17,