

## Datasheet

### Ly6g monoclonal antibody, clone 1A8 (PE)

**Catalog Number:** MAB15466

**Regulation Status:** For research use only (RUO)

**Product Description:** Rat monoclonal antibody raised against native mouse Ly6g.

**Clone Name:** 1A8

**Immunogen:** Ly6g-transfected EL4J cell line.

**Host:** Rat

**Theoretical MW (kDa):** 21-35

**Reactivity:** Mouse

**Applications:** Flow Cyt  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** PE

**Purification:** Affinity purification

**Isotype:** IgG2a

**Recommend Usage:** Flow Cytometry  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).

**Storage Instruction:** Store in the dark at 4°C. Avoid prolonged exposure to light.

**Entrez GeneID:** 546644

**Gene Symbol:** Ly6g

**Gene Alias:** Gr-1, Gr1, Ly-6G