

Datasheet

Ly6g monoclonal antibody, clone 1A8 (PerCP)

Catalog Number: MAB15467

Regulation Status: For research use only (RUO)

Product Description: Rat monoclonal antibody raised against native mouse Ly6g.

Clone Name: 1A8

Immunogen: Ly6g-transfected EL4J cell line.

Host: Rat

Theoretical MW (kDa): 21-35

Reactivity: Mouse

Applications: Flow Cyt
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PerCP

Purification: Affinity purification

Isotype: IgG2a

Recommend Usage: Flow Cytometry
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).

Storage Instruction: Store in the dark at 4°C. Avoid prolonged exposure to light.

Entrez GeneID: 546644

Gene Symbol: Ly6g

Gene Alias: Gr-1, Gr1, Ly-6G