

## Datasheet

### Fcgr3/Fcgr2b monoclonal antibody, clone 2.4G2 (FITC)

**Catalog Number:** MAB15473

**Regulation Status:** For research use only (RUO)

**Product Description:** Rat monoclonal antibody raised against native mouse Fcgr3/Fcgr2b.

**Clone Name:** 2.4G2

**Immunogen:** Mouse BALB/c macrophage J774 cell line.

**Host:** Rat

**Theoretical MW (kDa):** 40-60

**Reactivity:** Mouse

**Applications:** Flow Cyt

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** FITC

**Purification:** Affinity purification

**Isotype:** IgG2b

**Recommend Usage:** Flow Cytometry

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).

**Storage Instruction:** Store in the dark at 4°C. Avoid prolonged exposure to light.

**Entrez GeneID:** 14131|14130

**Gene Symbol:** Fcgr3

**Gene Alias:** CD16