

Datasheet

Fcgr3/Fcgr2b monoclonal antibody, clone 2.4G2 (PerCP)

Catalog Number: MAB15475

Regulation Status: For research use only (RUO)

Product Description: Rat monoclonal antibody raised against native mouse Fcgr3/Fcgr2b.

Clone Name: 2.4G2

Immunogen: Mouse BALB/c macrophage J774 cell line.

Host: Rat

Theoretical MW (kDa): 40-60

Reactivity: Mouse

Applications: Flow Cyt

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PerCP

Purification: Affinity purification

Isotype: IgG2b

Recommend Usage: Flow Cytometry

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).

Storage Instruction: Store in the dark at 4°C. Avoid prolonged exposure to light.

Entrez GeneID: 14131|14130

Gene Symbol: Fcgr3

Gene Alias: CD16