

Datasheet

Ptprc monoclonal antibody, clone 30-F11 (CF-Blue)

Catalog Number: MAB15485

Regulation Status: For research use only (RUO)

Product Description: Rat monoclonal antibody raised against native mouse Ptprc.

Clone Name: 30-F11

Immunogen: Mouse thymocytes and splenocytes.

Host: Rat

Theoretical MW (kDa): 180-240

Reactivity: Mouse

Applications: Flow Cyt
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: CF-Blue

Purification: Affinity purification

Isotype: IgG2b

Recommend Usage: Flow Cytometry
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).

Storage Instruction: Store in the dark at 4°C. Avoid prolonged exposure to light.

Entrez GeneID: 19264

Gene Symbol: Ptprc

Gene Alias: B220, CD45R, Cd45, Ly-5, Lyt-4, T200, loc