

## Datasheet

### Ptprc monoclonal antibody, clone 30-F11 (PerCP-Cyanine5.5)

**Catalog Number:** MAB15486

**Regulation Status:** For research use only (RUO)

**Product Description:** Rat monoclonal antibody raised against native mouse Ptprc.

**Clone Name:** 30-F11

**Immunogen:** Mouse thymocytes and splenocytes.

**Host:** Rat

**Theoretical MW (kDa):** 180-240

**Reactivity:** Mouse

**Applications:** Flow Cyt  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** PerCP-Cyanine5.5

**Purification:** Affinity purification

**Isotype:** IgG2b

**Recommend Usage:** Flow Cytometry  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).

**Storage Instruction:** Store in the dark at 4°C. Avoid prolonged exposure to light.

**Entrez GeneID:** 19264

**Gene Symbol:** Ptprc

**Gene Alias:** B220, CD45R, Cd45, Ly-5, Lyt-4, T200, loc