

Datasheet

Thy1.2 monoclonal antibody, clone 30-H12 (FITC)

Catalog Number: MAB15488

Regulation Status: For research use only (RUO)

Product Description: Rat monoclonal antibody raised against native mouse Thy1.2.

Clone Name: 30-H12

Immunogen: Mouse thymocytes and splenocytes.

Host: Rat

Theoretical MW (kDa): 25-35

Reactivity: Mouse

Applications: Flow Cyt

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: FITC

Purification: Affinity purification

Isotype: IgG2b

Recommend Usage: Flow Cytometry

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).

Storage Instruction: Store in the dark at 4°C. Avoid prolonged exposure to light.