

## Datasheet

### Thy1.2 monoclonal antibody, clone 30-H12 (PerCP)

**Catalog Number:** MAB15490

**Regulation Status:** For research use only (RUO)

**Product Description:** Rat monoclonal antibody raised against native mouse Thy1.2.

**Clone Name:** 30-H12

**Immunogen:** Mouse thymocytes and splenocytes.

**Host:** Rat

**Theoretical MW (kDa):** 25-35

**Reactivity:** Mouse

**Applications:** Flow Cyt

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** PerCP

**Purification:** Affinity purification

**Isotype:** IgG2b

**Recommend Usage:** Flow Cytometry

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).

**Storage Instruction:** Store in the dark at 4°C. Avoid prolonged exposure to light.