

## Datasheet

### Siglec1 monoclonal antibody, clone 3D6 (FITC)

**Catalog Number:** MAB15495

**Regulation Status:** For research use only (RUO)

**Product Description:** Rat monoclonal antibody raised against native mouse Siglec1.

**Clone Name:** 3D6

**Immunogen:** Native purified Siglec1 from spleen.

**Host:** Rat

**Theoretical MW (kDa):** 183

**Reactivity:** Mouse

**Applications:** Flow Cyt  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** FITC

**Purification:** Affinity purification

**Isotype:** IgG2a

**Recommend Usage:** Flow Cytometry  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).

**Storage Instruction:** Store in the dark at 4°C. Avoid prolonged exposure to light.

**Entrez GeneID:** 20612

**Gene Symbol:** Siglec1

**Gene Alias:** CD169, Siglec-1, Sn