

## Datasheet

### Emr1 monoclonal antibody, clone BM8 (PE)

TM7LN3

**Catalog Number:** MAB15509**Regulation Status:** For research use only (RUO)**Product Description:** Rat monoclonal antibody raised against native mouse Emr1.**Clone Name:** BM8**Immunogen:** Murine macrophages.**Host:** Rat**Theoretical MW (kDa):** 160**Reactivity:** Mouse**Applications:** Flow Cyt

(See our web site product page for detailed applications information)

**Protocols:** See our web site at<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols**Form:** Liquid**Conjugation:** PE**Purification:** Affinity purification**Isotype:** IgG2a**Recommend Usage:** Flow Cytometry

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).**Storage Instruction:** Store in the dark at 4°C. Avoid prolonged exposure to light.**Entrez GeneID:** 13733**Gene Symbol:** Emr1**Gene Alias:** DD7A5-7, EGF-TM7, F4/80, Gpf480, Ly71,