

Datasheet

Klrk1 monoclonal antibody, clone C7 (FITC)

Catalog Number: MAB15516

Regulation Status: For research use only (RUO)

Product Description: Hamster monoclonal antibody raised against native mouse Klrk1.

Clone Name: C7

Immunogen: Soluble mouse Klrk1.

Host: Hamster

Theoretical MW (kDa): 27

Reactivity: Mouse

Applications: Flow Cyt
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: FITC

Purification: Affinity purification

Isotype: IgG

Recommend Usage: Flow Cytometry
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).

Storage Instruction: Store in the dark at 4°C. Avoid prolonged exposure to light.

Entrez GeneID: 27007

Gene Symbol: Klrk1

Gene Alias: D6H12S2489E, NKG2-D, Nkg2d