

## Datasheet

### Klrk1 monoclonal antibody, clone C7 (PE)

**Catalog Number:** MAB15517

**Regulation Status:** For research use only (RUO)

**Product Description:** Hamster monoclonal antibody raised against native mouse Klrk1.

**Clone Name:** C7

**Immunogen:** Soluble mouse Klrk1.

**Host:** Hamster

**Theoretical MW (kDa):** 27

**Reactivity:** Mouse

**Applications:** Flow Cyt  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** PE

**Purification:** Affinity purification

**Isotype:** IgG

**Recommend Usage:** Flow Cytometry  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).

**Storage Instruction:** Store in the dark at 4°C. Avoid prolonged exposure to light.

**Entrez GeneID:** 27007

**Gene Symbol:** Klrk1

**Gene Alias:** D6H12S2489E, NKG2-D, Nkg2d