

## Datasheet

### Csf1r monoclonal antibody, clone CSF-1R (PE)

**Gene Alias:** AI323359, CD115, CSF-1R, Csfmr, Fim-2,  
Fms, M-CSFR

**Catalog Number:** MAB15521

**Regulation Status:** For research use only (RUO)

**Product Description:** Rat monoclonal antibody raised  
against recombinant mouse Csf1r.

**Clone Name:** CSF-1R

**Immunogen:** Recombinant protein corresponding to  
mouse Csf1r.

**Host:** Rat

**Theoretical MW (kDa):** 145

**Reactivity:** Mouse

**Applications:** Flow Cyt  
(See our web site product page for detailed applications  
information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product  
page for detailed protocols

**Form:** Liquid

**Conjugation:** PE

**Purification:** Affinity purification

**Isotype:** IgG2a

**Recommend Usage:** Flow Cytometry  
The optimal working dilution should be determined by  
the end user.

**Storage Buffer:** In PBS, pH 7.4 (protein stabilizer,  
0.09% sodium azide).

**Storage Instruction:** Store in the dark at 4°C. Avoid  
prolonged exposure to light.

**Entrez GeneID:** 12978

**Gene Symbol:** Csf1r