

Datasheet

Csf1r monoclonal antibody, clone CSF-1R (PerCP)

Gene Alias: AI323359, CD115, CSF-1R, Csfmr, Fim-2,
Fms, M-CSFR

Catalog Number: MAB15522

Regulation Status: For research use only (RUO)

Product Description: Rat monoclonal antibody raised
against recombinant mouse Csf1r.

Clone Name: CSF-1R

Immunogen: Recombinant protein corresponding to
mouse Csf1r.

Host: Rat

Theoretical MW (kDa): 145

Reactivity: Mouse

Applications: Flow Cyt
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Form: Liquid

Conjugation: PerCP

Purification: Affinity purification

Isotype: IgG2a

Recommend Usage: Flow Cytometry
The optimal working dilution should be determined by
the end user.

Storage Buffer: In PBS, pH 7.4 (protein stabilizer,
0.09% sodium azide).

Storage Instruction: Store in the dark at 4°C. Avoid
prolonged exposure to light.

Entrez GeneID: 12978

Gene Symbol: Csf1r