

## Datasheet

### Ly6a monoclonal antibody, clone D7 (PerCP)

TAP

**Catalog Number:** MAB15527**Regulation Status:** For research use only (RUO)**Product Description:** Rat monoclonal antibody raised against native mouse Ly6a.**Clone Name:** D7**Immunogen:** IL-2-dependent mouse T-cell line CTL-L.**Host:** Rat**Theoretical MW (kDa):** 18**Reactivity:** Mouse**Applications:** Flow Cyt

(See our web site product page for detailed applications information)

**Protocols:** See our web site at<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols**Form:** Liquid**Conjugation:** PerCP**Purification:** Affinity purification**Isotype:** IgG2a**Recommend Usage:** Flow Cytometry

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (protein stabilizer, 0.09% sodium azide).**Storage Instruction:** Store in the dark at 4°C. Avoid prolonged exposure to light.**Entrez GeneID:** 110454**Gene Symbol:** Ly6a**Gene Alias:** Ly-6A.2, Ly-6A/E, Ly-6E.1, Sca-1, Sca1,