

Datasheet

THY1 monoclonal antibody, clone CL1028

Catalog Number: MAB15660

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against partial recombinant human THY1.

Clone Name: CL1028

Immunogen: Recombinant protein corresponding to human THY1.

Sequence:

VTSLTACLVDQSLRLDCRHENTSSSPIQYEFSLTRETK
KHVLFGTVGVP EHTYRSRTNFTSKYNMKVLYLSAFTS
KDEGTYTCALHHS GHSPPISSQNVTVLRDKLVKCEGIS
LLAQNTSWL

Host: Mouse

Reactivity: Human

Applications: IHC-P, WB-Ti

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Protein A purification

Isotype: IgG2b

Recommend Usage: Immunohistochemistry

(Formalin/PFA-fixed paraffin-embedded sections)

(1:2500-1:5000)

Western Blot (1:500-1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (40% glycerol, 0.02% sodium azide).

Storage Instruction: Store at 4°C. For long term

storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 7070

Gene Symbol: THY1

Gene Alias: CD90, FLJ33325