

Datasheet

Cd34 monoclonal antibody, clone MEC14.7 (PURE)

Catalog Number: MAB15682

Regulatory Status: For research use only (RUO)

Product Description: Rat monoclonal antibody raised against mouse cd34.

Clone Name: MEC14.7

Immunogen: Mouse cd34

Host: Rat

Theoretical MW (kDa): 75-120

Reactivity: Mouse

Applications: Flow Cyt

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: PURE

Purification: Affinity chromatography

Isotype: IgG2a