

## Datasheet

### **Cd34 monoclonal antibody, clone MEC14.7 (Biotin)**

**Catalog Number:** MAB15685

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rat monoclonal antibody raised against mouse cd34.

**Clone Name:** MEC14.7

**Immunogen:** Mouse cd34

**Host:** Rat

**Theoretical MW (kDa):** 75-120

**Reactivity:** Mouse

**Applications:** Flow Cyt

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** Biotin

**Purification:** Affinity chromatography

**Isotype:** IgG2a