

Datasheet

p40 monoclonal antibody, clone L14-A (FITC)

Catalog Number: MAB16017

Regulation Status: For research use only (RUO)

Product Description: Rabbit monoclonal antibody raised against synthetic peptide of human p40.

Clone Name: L14-A

Host: Rabbit

Reactivity: Human

Applications: Flow Cyt
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: FITC

Purification: EVAC purification

Isotype: IgG

Recommend Usage: Flow Cytometry
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (10 mg/mL BSA, 0.05% Sodium Azide).

Storage Instruction: Store in the dark at 4°C. Avoid prolonged exposure to light. Do not freeze.