

Datasheet

Cd69 monoclonal antibody, clone H1.2F3 (FITC)

Gene Symbol: Cd69

Gene Alias: 5830438K24Rik, AI452015, AIM, VEA

Catalog Number: MAB16150

Regulatory Status: For research use only (RUO)

Product Description: Hamster monoclonal antibody raised against mouse Cd69.

Clone Name: H1.2F3

Immunogen: Y245 mouse dendritic epidermal T cell line.

Host: Hamster

Reactivity: Mouse

Applications: Flow Cyt, IHC-P, IP
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: FITC

Purification: Purified

Isotype: IgG, lamda

Recommend Usage: Flow Cytometry
Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) (10 ug/mL)
Immunoprecipitation
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.09% sodium azide).

Storage Instruction: Store at 4°C. For long term storage store at -20°C.
Avoid prolonged exposure to light and repeated freezing and thawing.

Entrez GeneID: 12515