

Datasheet

Clcn3 monoclonal antibody, clone S258-5 (FITC)

Gene Symbol: Clcn3

Gene Alias: CIC-3

Catalog Number: MAB16803

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against synthetic peptide of rat Clcn3.

Clone Name: S258-5

Immunogen: A synthetic peptide corresponding to amino acids 98-115 at N-terminus of rat Clcn3.

Host: Mouse

Reactivity: Human, Rat

Applications: ICC, IF, IHC-P, WB-Ti
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: FITC

Purification: Protein G purification

Isotype: IgG1

Recommend Usage: Immunocytochemistry (1:100)
Immunofluorescence (1:100)
Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) (1:1000)
Western Blot (1:1000)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide).

Storage Instruction: Store at -20°C.

Entrez GeneID: 84360