

Datasheet

Git1 monoclonal antibody, clone S39B-8 (FITC)

Catalog Number: MAB16812

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against partial recombinant rat Git1.

Clone Name: S39B-8

Immunogen: Recombinant protein corresponding to amino acids 375-770 at C-terminus of rat Git1.

Host: Mouse

Reactivity: Human, Rat

Applications: ICC, IF, IP, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: FITC

Purification: Protein G purification

Isotype: IgG1

Recommend Usage: Immunocytochemistry

Immunofluorescence

Immunoprecipitation

Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide).

Storage Instruction: Store at -20°C.

Entrez GeneID: 83709

Gene Symbol: Git1

Gene Alias: -

References:

1. The multifunctional GIT family of proteins. Hoefen RJ, Berk BC. J Cell Sci. 2006 Apr 15;119(Pt 8):1469-75.