

Datasheet

Accn2 monoclonal antibody, clone S271-44 (FITC)

Catalog Number: MAB16836

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against partial recombinant mouse Accn2.

Clone Name: S271-44

Immunogen: Recombinant protein corresponding to amino acids 460-526 at C-terminus of mouse Accn2.

Host: Mouse

Reactivity: Human, Mouse, Rat

Applications: ICC, IF, IHC-P, WB-Ti
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: FITC

Purification: Protein G purification

Isotype: IgG1

Recommend Usage: Immunocytochemistry
Immunofluorescence
Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections)
Western Blot (1:1000)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide).

Storage Instruction: Store at -20°C.

Entrez GeneID: 11419

Gene Symbol: Accn2

Gene Alias: AI843610, ASIC, ASIC1, ASIC1a, B530003N02Rik, BNaC2