

Datasheet

Bcan monoclonal antibody, clone S294A-6 (FITC)

Gene Alias: ALPBRE

Catalog Number: MAB16839

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against partial recombinant rat Bcan.

Clone Name: S294A-6

Immunogen: Recombinant protein corresponding to amino acids 219-655 of rat Bcan.

Host: Mouse

Reactivity: Human, Rat

Applications: ICC, IF, WB-Ti

(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: FITC

Purification: Protein G purification

Isotype: IgG2b

Recommend Usage: Immunocytochemistry

Immunofluorescence

Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide).

Storage Instruction: Store at -20°C.

Entrez GeneID: 25393

Gene Symbol: Bcan