

Datasheet

Mfn2 monoclonal antibody, clone S153-5 (FITC)

Gene Symbol: Mfn2

Gene Alias: D630023P19Rik, Fzo

Catalog Number: MAB16872

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against partial recombinant human Mfn2.

Clone Name: S153-5

Immunogen: Recombinant protein corresponding to amino acids 370-600 at N-terminus of mouse Mfn2.

Host: Mouse

Reactivity: Human, Mouse, Rat

Applications: ICC, IF, IHC-P, WB-Ti

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: FITC

Purification: Protein G purification

Isotype: IgG2a

Recommend Usage: Immunocytochemistry (1:100)

Immunofluorescence (1:100)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections)

Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide).

Storage Instruction: Store at -20°C.

Entrez GeneID: 170731