

Datasheet

Nlgn3 monoclonal antibody, clone S110-29 (FITC)

Gene Symbol: Nlgn3

Gene Alias: -

Catalog Number: MAB16893

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against partial recombinant rat Nlgn3.

Clone Name: S110-29

Immunogen: Recombinant protein corresponding to amino acids 730-848 at C-terminus of rat Nlgn3.

Host: Mouse

Reactivity: Human, Mouse, Rat

Applications: ICC, IF, IHC-P, WB-Ti
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: FITC

Purification: Protein G purification

Isotype: IgG1

Recommend Usage: Immunocytochemistry
Immunofluorescence
Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections)
Western Blot (1:1000)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide).

Storage Instruction: Store at -20°C.

Entrez GeneID: 171297