

Datasheet

4-Hydroxynonenal monoclonal antibody, clone 12F7 (FITC)

Catalog Number: MAB17401

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against synthetic 4-Hydroxynonenal (4-HNE).

Clone Name: 12F7

Immunogen: Synthetic 4-Hydroxynonenal modified Keyhole Limpet Hemocyanin (KLH).

Host: Mouse

Applications: ELISA, ICC, IF, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: FITC

Purification: Protein G purification

Isotype: IgG1

Recommend Usage: ELISA (1:1000)

Immunocytochemistry (1:50)

Immunofluorescence (1:50)

Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide).

Storage Instruction: Store at -20°C.