

## Datasheet

### Malondialdehyde monoclonal antibody, clone 6H6 (FITC)

**Catalog Number:** MAB17409

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against synthetic Malondialdehyde (MDA).

**Clone Name:** 6H6

**Immunogen:** Synthetic Malondialdehyde modified Keyhole Limpet Hemocyanin (KLH).

**Host:** Mouse

**Applications:** ELISA, Flow Cyt, ICC, IF, WB  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** FITC

**Purification:** Protein G purification

**Isotype:** IgG1

**Recommend Usage:** ELISA (1:1000)

Flow Cytometry (1:50)

Immunocytochemistry (1:50)

Immunofluorescence (1:50)

Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide).

**Storage Instruction:** Store at -20°C.